



ANTI Xa ASSAY (FOR MONITORING LOW MOLECULAR WEIGHT HEPARIN)

Synonym(s):

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL blood containing 0.3 mL of 3.2% sodium citrate
Method	Chromogenic
Reference Interval/Value	Therapeutic range for Enoxaparin (once daily dosing): 1.0 - 2.0 IU/ml Enoxaparin (twice daily dosing): 0.5 - 1.0 IU/ml Tinzaparin: 0.85 - 1.0 IU/ml
Turnaround Time	3 – 4 days
Day(s) Test Set up	Monday – Saturday
Remarks	Please indicate the type of heparin given.

Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608