

Make/Change Appointment Pay Online

## CREATINE KINASE (CK), SERUM

Synonym(s): CK

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Clinical Biochemistry  Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	3 mL plain blood. Avoid haemolysis.  Refrain from exercise or strenuous activity prior to blood collection.
Method	Enzymatic
Reference Interval/Value	1 - 30 days : 5 - 183 U/L 31 - 365 days : 5 - 146 U/L 1 - 6 yrs : 5 - 163 U/L 7 - 9 yrs : 5 - 177 U/L 10 - 12 yrs : 5 - 217 U/L 13 - 17 yrs : 5 - 251 U/L ≥ 18 yrs : M 56 - 336; F 44 - 201 U/L
Turnaround Time	Mean = 4 hrs, Range: 3 - 5 hrs
Day(s) Test Set up	Daily, 24 hrs STAT

Back To Previous Page (javascript: history.go(-1))

## **Change History Notes**

29 Apr 2013 06:37 PM
 Reference Interval has been last updated on 22 Apr 2013.

https://www.sgh.com.sg/Clinical-Departments-Centers/Pathology/Pathology/Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/... 19/10/2018

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608