



## **ELECTRON MICROSCOPY**

**Synonym(s):** Ultrastructural Studies

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Histopathology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/histopathology.aspx).
Description	Specimen is placed in 2.5% glutaraldehyde in cacodylate buffer (the time for which the tissue is placed into fixative must be recorded).  Contact the laboratory (Tel: 6321 4008) for further information, if necessary.

## Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608