



GLUCOSE, SERUM

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	3 mL plain blood. Separate serum within 2h of collection. Send specimen immediately to lab. Collect in sodium fluoride tube if sample will take more than 2h to reach the lab.
Method	Enzymatic/Oxygen Electrode
Reference Interval/Value	From 0 day Fasting : 3.9 - 6.0 mmol/L ; Random : 3.9 - 11.0 mmol/L
Turnaround Time	Mean = 4 hrs, Range: 3 – 5 hrs
Day(s) Test Set up	Daily, 24 hrs STAT

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

