>



 Make/Change Appointment
 S
 Pay Online

KAPPA, LAMBDA FREE LIGHT CHAINS, SERUM

Synonym(s):

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology- handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	3 mL plain blood
Method	Immunoassay
Reference Interval/Value	From 0 day: kappa : 3.3 – 19.4 mg/L lambda : 5.7 – 26.3 mg/L kappa/lambda ratio : 0.26 – 1.65
Turnaround Time	Mean = 7 days, Range: 3 - 11 days
Day(s) Test Set up	Thursday

Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

<