

V₈ \$ Make/Change Appointment

Pay Online

METHAEMOGLOBIN (QUANTITATIVE) BLOOD

Synonym(s):

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments- centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical- biochemistry.aspx).
Specimen Required	3 mL blood in heparinised syringe immersed in ice and sent within 2 h after collection. Sodium fluoride/potassium oxalate tubes are not suitable for analysis.
Method	Co-oximetry
Reference Interval/Value	from 0 day : 0 – 1.5%
Turnaround Time	Mean = 4 hrs, Range: 3 – 5hrs
Day(s) Test Set up	Daily, 24 hrs STAT

Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608 <