



## PROTEIN C ASSAY (FUNCTIONAL)

## Synonym(s):

<

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL blood containing 0.3 mL of 3.2% sodium citrate
Method	Chromogenic Assay
Reference Interval/Value	70 - 140%
Turnaround Time	1 week
Day(s) Test Set up	Weekly
Remarks	<b>Note</b> : Coumadin will lower the Protein C level. Assay cannot be performed if patient is on heparin.

## Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608