



pH, URINE

Synonym(s):

<

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments- centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical- biochemistry.aspx).
Specimen Required	20 mL fresh random urine in sterile container.
Method	H+ Ion-selective Electrode Potentiometry
Reference Interval/Value	4.6 - 7.0
Turnaround Time	Mean = 4 hrs, Range: 3 - 5 hrs
Day(s) Test Set up	Monday - Saturday

Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608