



## ACANTHAMOEBA/NAEGLERIA CULTURE

**Synonym(s):** Naegleria/Acanthamoeba Culture

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Bacteriology/Mycology/Parasitology - Parasites & Other Investigations Click here to find out more about the write-up (/Clinical-Departments- Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/pages/bacteriology-mycology-parasitology.aspx).
Specimen Required	Corneal scraping, cerebrospinal fluid, brain biopsy tissue
Method	Culture on non-nutrient agar containing Escherichia coli
Test Result	Positive for Acanthamoeba/Naegleria fowleri; No Acanthamoeba/Naegleria fowleri isolated
Turnaround Time	10 days
Day(s) Test Set up	Monday - Saturday

## Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608