



ACTIVATED PARTIAL THROMBOPLASTIN TIME (APTT)

Synonym(s): APTT, PTT

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL blood containing 0.3 mL of 3.2% sodium citrate
Method	Photo-optical Clot Detection
Reference Interval/Value	25.7 – 32.9 seconds (CS 2100)
Turnaround Time	3 – 4 hours
Day(s) Test Set up	Daily

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Change History Notes

- **14 May 2014 02:33 PM**
Reference Interval has been last updated on 14 May 2014.