



	Make/Change Appointment		Pay Online
---	---	---	----------------------------

AMINO ACIDS, CSF - test referred to KKH Biochemical Genetics Lab w.e.f 03 October 2016

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx) .
Specimen Required	200 uL CSF collected in a STERILE CONTAINER (this test should be ordered in conjunction with Amino Acids, Plasma (lithium heparin tube). The specimens for both tests (Amino Acids, Plasma and CSF) should be collected within 1 hour and pack in ice pack. (KKH requirement)
Storage and Transport	Keep cold with ice pack and send immediately.

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

