



ANTI-DOUBLE STRAND DNA IIF (ANTI-CRITHIDIA ANTIBODY)

Synonym(s): Anti DS DNA (IF), ds DNA (IIF) Antibody, Anti Crithidia Luciliae Antibody, Crithidia Antibody

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Immunology/Serology/STD/Allergy - Immunology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/immunologyserologystdallergy.aspx).
Specimen Required	3 mL plain blood
Method	Immunofluorescence
Test Result	Screening dilution is 1:10. Reported as Positive or Negative
Turnaround Time	Range 1 – 2 days
Day(s) Test Set up	Monday - Friday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

