



	Make/Change Appointment		Pay Online
--	-------------------------	--	------------

ANTI-EXTRACTABLE NUCLEAR ANTIGENS (ENA)-PROFILE

Synonym(s): Anti ENA, ENA Profile, Anti SSA, Anti Ro, Anti SSB, Anti La, Anti Sm, Anti Smith, Anti RNP, Anti Jo-1, Anti Scl-70

Back To Previous Page (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>) (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>)

Lab Section Category	Immunology/Serology/STD/Allergy - Immunology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/immunologyserologystdallergy.aspx).
Specimen Required	3 mL plain blood
Method	Enzyme Immunoassay
Test Result	Reported as numerical readings of ratio with interpretations for the auto-antibodies to: - Anti-Sm - Anti-RNP - Anti-SSA (Ro) - Anti-SSB (La) - Anti-Scl-70 - Anti-Jo-1
Turnaround Time	Range 2 – 4 days
Day(s) Test Set up	Monday to Friday

Back To Previous Page (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>) (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>)

Change History Notes

- **05 Apr 2018 03:03 PM**
Method has been updated w.e.f. 05 March 2018