



ANTI-GAD ANTIBODY ELISA

Synonym(s): GAD, Glutamic Acid Decarboxylase Antibody

Back To Previous Page (<http://cms-sgh.sppub/clinical-departments-centers/pathology/pathology-handbook/pages/testlisting.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAtIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>) (<http://cms-sgh.sppub/clinical-departments-centers/pathology/pathology-handbook/pages/testlisting.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAtIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>)

Lab Section Category	Immunology/Serology/STD/Allergy - Immunology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/immunologyserologystdallergy.aspx).
Specimen Required	3 mL plain blood
Method	Enzyme Immunoassay
Test Result	Reported as numerical readings in IU/mL with interpretation.
Turnaround Time	Range 2-8 days
Day(s) Test Set up	Batch tested, once/week

Back To Previous Page (<http://cms-sgh.sppub/clinical-departments-centers/pathology/pathology-handbook/pages/testlisting.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAtIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>) (<http://cms-sgh.sppub/clinical-departments-centers/pathology/pathology-handbook/pages/testlisting.aspx?category=SW1tdW5vbG9neS9TZXJvbG9neS9TVEQvQWxsZXJneSAtIEltbXVub2xvZ3k%3d-qLAjr1Ja5yM%3d>)

Change History Notes

- **12 Mar 2019 12:19 PM**
New test available w.e.f 04 March 2019