



ANTI-SMOOTH MUSCLE ANTIBODY

Synonym(s): Anti Smooth Muscle, Smooth Muscle Antibody, Smooth Muscle Ab, SMA

Back To Previous Page ([javascript: history.go\(-1\)](#))

Lab Section Category	Immunology/Serology/STD/Allergy - Immunology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/immunologyserologystdallergy.aspx).
Specimen Required	3 mL plain blood
Method	Immunofluorescence
Test Result	Starting dilution is 1:100. Serum worked-up in doubling dilutions from 1:100 to 1:400. Results reported as Negative, or Positive with titre.
Turnaround Time	Range 2 – 4 days
Day(s) Test Set up	Monday – Friday

Back To Previous Page ([javascript: history.go\(-1\)](#))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

