



ANTI-Xa (RIVAROXABAN)

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL blood containing 0.3mL of 3.2% sodium citrate
Method	Chromogenic
Reference Interval/Value	There are no established therapeutic ranges for rivaroxaban. Please note that plasma rivaroxaban levels are affected by the timing of blood sampling and renal function.
Turnaround Time	3 – 4 days
Day(s) Test Set up	Monday – Saturday (office hours)
Remarks	Nil

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

