

Make/Change Appointment Pay Online

AVIAN INFLUENZA H5N1 PCR (QUALITATIVE)

Synonym(s): H5, H5N1

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Molecular Diagnostics - Microbiology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/moleculardiagnostics.aspx).
Indications	For detection of the avian influenza A/H5N1 in patients with severe respiratory infection and travel to or residence in a country where avian influenza A/H5N1 has been detected. The influenza multiplex PCR test must be ordered with this test.
Specimen Required	 a) Throat/nasopharyngeal swabs (preferred sample for adults): Use only Dacron-tipped or flocked swab. Send dry in sterile container or swab sheath. b) Nasopharyngeal aspirate/wash, endotracheal aspirate, bronchoalveolar lavage: Send at least 1 mL of sample in sterile leak-proof screw-capped container without any media. Paste red bioharzard alert stickers on clinical sample bags for SGH inpatients.
Storage and Transport	Refrigerate sample until transfer to laboratory. Do not freeze . Send sample at 2 - 8°C (with an ice-pack). Sample must reach laboratory within 24 hours after collection. Please contact Molecular Laboratory (63266920) before sending sample.
Method	Reverse-transcription Polymerase Chain Reaction (In-house Real-time)
Test Result	Detected, Not Detected or Inconclusive

Turnaround Time	Urgent SGH cases only: 6 hours (requires approval of microbiologist) Routine cases: 1 working day Patient samples for urgent/stat testing must be from patients who meet MOH criteria of a suspect case. These must be ordered by an infectious diseases physician and discussed with a Molecular Laboratory microbiologist.
Day(s) Test Set up	Monday-Saturday
Remarks	In addition to universal precautions, strict barrier precautions and <u>respiratory protection</u> are required in the collection of respiratory samples from patients suspected to be infected with avian influenza. This includes the use of personal protective equipment like gloves, gowns and apron, safety eyewear and powered air purifying respirators (PAPR) or N95 masks.

Back To Previous Page (javascript: history.go(-1))

Change History Notes

- 19 May 2014 04:11 PM
 Specimen Required: Specimen Required has been last updated on 19/05/14.
- 04 Dec 2015 12:16 PM
 Storage and Transportation: Storage and Transportation has been last updated on 04/12/15.

Turnaround Time: Turnaround Time has been updated on 04/12/15.

• 16 Jan 2017 10:11 AM

<

Specimen Required: Specimen Required has been last updated on 16/01/17

Turnaround Time: Turnaround Time has been last updated on 16/01/17

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608