



## BETA-2-MICROGLOBULIN/CREATININE PROFILE, URINE

### Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

<b>Lab Section Category</b>	<b>Clinical Biochemistry</b> Click here to find out more about the <b>write-up</b> ( <a href="#">/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx</a> ).
<b>Specimen Required</b>	20 mL urine in sterile container  1. Void first morning urine and discard. 2. Drink one glass of water. 3. Collect urine 1h later. Send 20 mL immediately to lab
<b>Method</b>	Immunoassay /Spectrophotometry
<b>Reference Interval/Value</b>	26 - 260 µg/g creatinine
<b>Turnaround Time</b>	Mean = 4 hrs, Range: 3 - 5 hrs
<b>Day(s) Test Set up</b>	Monday - Saturday
<b>Remarks</b>	Measured Parameters: Beta-2-microglobulin (Urine) and Creatinine (Urine)

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

