



BLOOD FOR MICROFILARIA

Synonym(s): Microfilaria, Blood

Back To Previous Page ([javascript: history.go\(-1\)](#))

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL EDTA blood or thick films To be taken preferably around midnight.
Method	Microscopy
Test Result	Absent
Turnaround Time	3 – 4 hours
Day(s) Test Set up	Monday – Saturday (office hours)

Back To Previous Page ([javascript: history.go\(-1\)](#))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

