



## BONE MARROW TREPINE BIOPSIES

### Synonym(s):

**Back To Previous Page** (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SGlzdG9wYXRob2xvZ3k%3d-%2f%2fUm1oSD06M%3d>) (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SGlzdG9wYXRob2xvZ3k%3d-%2f%2fUm1oSD06M%3d>)

<b>Lab Section Category</b>	<b>Histopathology</b> Click here to find out more about the <b>write-up</b> ( <a href="/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/histopathology.aspx">/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/histopathology.aspx</a> ).
<b>Description</b>	<ol style="list-style-type: none"> <li>1. Specimens should be placed in 10% buffered formalin.</li> <li>2. Label containers and complete Histo CPOE order.</li> </ol>

**Back To Previous Page** (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SGlzdG9wYXRob2xvZ3k%3d-%2f%2fUm1oSD06M%3d>) (<http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=SGlzdG9wYXRob2xvZ3k%3d-%2f%2fUm1oSD06M%3d>)

©Copyright 2019 Singapore General Hospital, Outram Road, Singapore 169608

