



	Make/Change Appointment		Pay Online
---	---	---	----------------------------

CADMIUM, BLOOD

Synonym(s):

[Back To Previous Page](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up .
Specimen Required	5 mL blood collected in K2EDTA trace-metal free tube. Do not centrifuge as whole blood is required.
Method	Inductively Coupled Plasma Mass Spectrometry
Reference Interval/Value	≥ 20 yrs: <2.7 ug/L
Turnaround Time	Mean = 7 days, Range 3 - 11 days
Day(s) Test Set up	Monday - Friday
Remarks	Instructions Before Sampling: Patient should avoid taking either gadolinium or iodine containing contrast media for at least 4 days as they can interfere with trace metal tests.

[Back To Previous Page](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

