



CD 4/8 ASSAY

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).			
Specimen Required	3 mL EDTA blood			
Method	Flow Cytometry			
Reference Interval/Value	Cell markers	Unit	≤ 50 years	> 50 years
	CD 4	%	28.2 – 50.7	29.8 – 52.3
	CD 4 absolute #	X10 ⁹ /L	288 – 1507	256 – 1594
	CD 8	%	12.5 – 38.5	8.8 – 38.6
	CD 8 absolute #	X10 ⁹ /L	148 – 1039	59 – 1034
	CD 4/CD 8 ratio		0.7 – 2.5	0.7 – 2.5
Turnaround Time	1 – 2 days			
Day(s) Test Set up	Monday – Saturday (office hours)			

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

