



	<a href="#">Make/Change Appointment</a>		<a href="#">Pay Online</a>
---	---	---	----------------------------

## CHIKUNGUNYA/DENGUE DUAL ISOLATION

**Synonym(s):** Chikungunya / Dengue Duo Culture

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

<b>Lab Section Category</b>	<b>Virology</b> Click here to find out more about the <b>write-up</b> ( <a href="/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx">/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx</a> ).
<b>Specimen Required</b>	6 mL EDTA blood or 3-5 mL plain blood during the first week of illness.  (See also <b>Virology Section, General Guidelines for Virus Isolation &amp; Collection Methods for Virus Isolation by Sample Type</b> ( <a href="http://shintersp.shses.shs.com.sg:82/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx">http://shintersp.shses.shs.com.sg:82/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx</a> ) ( <a href="http://shintersp.shses.shs.com.sg:82/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx">http://shintersp.shses.shs.com.sg:82/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx</a> ) )
<b>Storage and Transport</b>	Refrigerate samples until transport to laboratory. <b>Do not freeze.</b>  Send samples on ice pack to laboratory.
<b>Method</b>	Tissue culture
<b>Test Result</b>	Reported as Positive or Negative
<b>Turnaround Time</b>	5 - 8 days
<b>Day(s) Test Set up</b>	Monday - Saturday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

### Change History Notes

- **02 May 2018 12:54 PM**  
Changed test name and TAT.

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608