



## D-DIMER TEST QUANTITATIVE

### Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

<b>Lab Section Category</b>	<b>Haematology</b> Click here to find out more about the <b>write-up</b> ( <a href="#">/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx</a> ).
<b>Specimen Required</b>	3 mL blood containing 0.3 mL of 3.2% sodium citrate
<b>Method</b>	Immuno-turbidimetry
<b>Reference Interval/Value</b>	0.19 – 0.55 mg/L FEU
<b>Turnaround Time</b>	1 day
<b>Day(s) Test Set up</b>	Monday - Sunday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

