



DENGUE VIRUS ANTIGEN

Synonym(s): Dengue Ag

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Virology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx).
Specimen Required	3 – 5 mL plain blood, within 9 days of onset of illness
Storage and Transport	Refrigerate samples if delay in transporting to laboratory is anticipated. Do not freeze.
Method	Enzyme immunoassay
Test Result	Reported as Positive or Negative
Turnaround Time	<1 – 3 days
Day(s) Test Set up	Monday, Wednesday and Friday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Change History Notes

- **29 Jul 2014 04:43 PM**
Change of Day(s) Test Setup to every Monday, Wednesday and Friday w.e.f 1 August 2014.