



EOSINOPHIL COUNT

Synonym(s):

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments- centers/pathology/pathology-handbook/lab-discipline-special- instructions/pages/haematology.aspx).
Specimen Required	3 mL EDTA blood
Method	Microscopy
Reference Interval/Value	$(0.04 - 0.44) \times 10^9/L$
Turnaround Time	4 hours
Day(s) Test Set up	Monday - Saturday (office hours)

Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608