



## **FLT 3 ITD MUTATIONS**

## Synonym(s):

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments- centers/pathology/pathology-handbook/lab-discipline-special- instructions/pages/haematology.aspx).
Specimen Required	10 mL EDTA bone marrow or 15 mL EDTA blood
Method	Polymerase Chain Reaction
Test Result	Absent
Turnaround Time	4 weeks
Day(s) Test Set up	Monday - Saturday (office hours)

## Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608