



HAEMOGLOBIN, UNSTABLE

Synonym(s):

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL EDTA blood
Method	Isopropanol Solubility
Test Result	Absent
Turnaround Time	1 day
Day(s) Test Set up	Monday - Friday (office hours)
Remarks	By appointment only. Note: False positive results will be obtained in blood specimens that contain > 5% fetal haemoglobin or in specimens that are received more than a week after the blood has been drawn.

Back To Previous Page (javascript: history.go(-1))

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608