



HAEMOPHILUS DUCREYI (CHANCROID) CULTURE

Synonym(s): Chancroid (Haemophilus Ducreyi) Culture

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Bacteriology/Mycology/Parasitology - Bacteria Click here to find out more about the write-up (/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/pages/bacteriology-mycology-parasitology.aspx) .
Specimen Required	Swab from genital ulcers, bubo aspirates and cervical swabs must be inoculated directly onto culture plates in the form of a "Z" at the bedside. Culture plates can be obtained from Client Services (Tel : 6321 4950/4952/4904) on request.
Storage and Transport	The plates must be despatched to the laboratory as soon as possible.
Method	Culture
Test Result	<i>Haemophilus ducreyi</i> reported No <i>Haemophilus ducreyi</i> isolated
Turnaround Time	Range 5 – 7 days
Day(s) Test Set up	Monday – Saturday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

