



HISTOPLASMA CAPSULATUM ANTIBODY TEST

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Bacteriology/Mycology/Parasitology - Fungi Click here to find out more about the write-up (/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/pages/bacteriology-mycology-parasitology.aspx).
Specimen Required	3 mL blood in plain tube / 0.5 mL cerebrospinal fluid
Method	Micro-immunodiffusion
Test Result	Positive for <i>Histoplasma</i> Ab M and H band; Positive for <i>Histoplasma</i> Ab M band; Negative for <i>Histoplasma</i> Ab M and H bands
Turnaround Time	24 hours
Day(s) Test Set up	Once a week

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

