



KAPPA, LAMBDA FREE LIGHT CHAINS, URINE

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	20 mL fresh urine in sterile container.
Method	Immunoassay
Reference Interval/Value	From 0 day: kappa : 0.78 - 13.5 mg/L lambda : 2.2 – 5.9 mg/L kappa/lambda ratio : 0.46 – 4.00
Turnaround Time	Mean = 7 days, Range: 3 – 11 days
Day(s) Test Set up	Thursday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

