



LACTATE, CSF

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	1 mL CSF in sodium fluoride/potassium oxalate as preservative. Keep cold with ice and send immediately.
Method	Spectrophotometry
Reference Interval/Value	Not available Contact Clinical Biochemistry Laboratory for assistance where necessary.
Turnaround Time	Mean = 4 hrs, Range: 3 – 5 hrs
Day(s) Test Set up	Daily

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Change History Notes



- **10 Feb 2017 12:01 PM**

Reference Interval has been last updated on 10 Feb 2017.

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

