



LACTATE, PLASMA

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	3 mL sodium fluoride/potassium oxalate blood. Keep cold with ice and send immediately.
Method	Spectrophotometry
Reference Interval/Value	From 0 day : 0.5 - 2.2 mmol/L
Turnaround Time	Mean = 4 hrs, Range: 3 - 5 hrs
Day(s) Test Set up	Daily

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

