



LEAD, BLOOD

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	5 mL blood collected in K2EDTA trace-metal free tube Do not centrifuge as whole blood is required.
Method	Inductively Coupled Plasma Mass Spectrometry
Reference Interval/Value	From 0 day: <100ug/L
Turnaround Time	Mean = 7 days ; Range 3 - 11 days
Day(s) Test Set up	Monday - Friday

Remarks	Instructions Before Sampling: <ol style="list-style-type: none">1. Patient should avoid taking either gadolinium or iodine containing contrast media for at least 4 days as they can interfere with trace metal tests.2. Finger prick method is discouraged for sampling as it can contribute to contamination.
----------------	---

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608