



## NICKEL, URINE

### Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

<b>Lab Section Category</b>	<b>Clinical Biochemistry</b> Click here to find out more about the <b>write-up</b> ( <a href="/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx">/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx</a> ).
<b>Specimen Required</b>	20 mL fresh random urine in sterile container
<b>Method</b>	Inductively Coupled Plasma Mass Spectrometry
<b>Reference Interval/Value</b>	≥ 20 yrs: <14 ug/L
<b>Turnaround Time</b>	Mean = 7 days ; Range 3 - 11 days
<b>Day(s) Test Set up</b>	Monday - Friday
<b>Remarks</b>	Instructions Before Sampling: Patient should avoid taking either gadolinium or iodine containing contrast media for at least 4 days as they can interfere with trace metal tests.

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

