



	Make/Change Appointment		Pay Online
---	---	---	----------------------------

NOROVIRUS ANTIGEN

Synonym(s): Norovirus Ag

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Virology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/virology.aspx) .
Specimen Required	1 ml of liquid stool or cherry-size solid stool in a clean plain screw-capped container, collected preferably within 3 days of onset of illness.
Storage and Transport	Refrigerate samples until transport to laboratory. Do not freeze. Transport samples on ice pack to laboratory.
Method	Enzyme Immunoassay
Test Result	Reported as Negative or Positive
Turnaround Time	<1 -7 days
Day(s) Test Set up	Tuesday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

