



OLIGOCLONAL BANDS, CSF

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	3 mL CSF in sterile container and 3mL plain blood
Method	Agarose Gel Isoelectric Focusing
Reference Interval/Value	Not detected
Turnaround Time	Mean = 12 days, Range: 8 – 19 days
Day(s) Test Set up	Monday – Friday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

