



OSMOLALITY, SERUM

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Clinical Biochemistry Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	3 mL plain blood
Method	Freezing Point Depression Osmometry
Reference Interval/Value	From 0 day : 275 - 301 mmol/kg
Turnaround Time	Mean = 4 hrs, Range: 3 – 5 hrs
Day(s) Test Set up	Daily

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Change History Notes

- **12 Mar 2012 06:38 PM**
Method has been last updated on 28 Feb 12.

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

