



PLATELET FUNCTION TEST

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	6 tubes of 3 mL blood containing 0.3 mL of 3.2% sodium citrate.
Storage and Transport	Specimen should be kept at room temperature and should not be kept in ice.
Reference Interval/Value	Platelet Aggregation with: (ADP) 64 – 111% (Collagen) 68 – 117% (Epinephrine) 46 – 122% (Ristocetin) 80 – 115% (Arachidonic Acid) 52 – 110%
Turnaround Time	1 day
Day(s) Test Set up	Monday - Friday

Remarks	Includes aggregation with 5 reagents. Note: To make appointment with laboratory (Tel: 6326 6022). Patient to fast overnight prior to test and should not be on aspirin for at least 2 weeks
----------------	--

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Change History Notes

- **14 May 2014 02:59 PM**
The Remarks section has been last updated on 14 May 2014.

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

