



	Make/Change Appointment		Pay Online
---	---	---	----------------------------

PROTHROMBIN TIME

Synonym(s): PT

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL blood containing 0.3 mL of 3.2% sodium citrate
Method	Photo-optical Clot Detection
Reference Interval/Value	9.9 – 11.4 seconds (CS2100)
Turnaround Time	3 - 4 hours
Day(s) Test Set up	Daily
Remarks	INR will be given upon request

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Change History Notes

- 14 May 2014 03:19 PM**
 Reference Interval has been last updated on 14 May 2014.

