



RETICULOCYTE COUNT

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	3 mL EDTA blood
Method	Supravital Staining with Brilliant Cresyl Blue
Reference Interval/Value	1 – 8 days : 2.5 – 6.5% > 8 days : 0.2 – 2.0%
Turnaround Time	2 - 3 hours
Day(s) Test Set up	Monday - Sunday

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

