



SOLUBLE FIBRIN/FIBRINOGEN MONOMER

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx) .
Specimen Required	3 mL blood containing 0.3 mL of 3.2% sodium citrate
Method	Haemagglutination
Test Result	Absent
Turnaround Time	3 - 4 hours
Day(s) Test Set up	Monday - Saturday (office hours)

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

