



URINE HAEMOSIDERIN

Synonym(s):

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

Lab Section Category	Haematology Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/haematology.aspx).
Specimen Required	12 mL fresh urine
Method	Prussian Blue Reaction
Test Result	Negative
Turnaround Time	4 hours
Day(s) Test Set up	Monday - Saturday (office hours)

[Back To Previous Page \(javascript: history.go\(-1\)\)](#)

©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608

