SINGAPORE GENERAL HOSPITAL DEPARTMENT OF CLINICAL PATHOLOGY CLIENT & SPECIMEN MANAGEMENT

REQUEST FOR SPECIMEN CONTAINERS AND COLLECTION MEDIA

(For SGH Campus Clients)

Requesting Staff:
Phone number :

Quantity To Be Entered By Requestor				
Item Required	Quantity	Item Required	Quantity	
BACTERIOLOGY SECTION		VIROLOGY SECTION		
Blood Culture Bottle (Aerobic)		Hank's Transport Medium		
Blood Culture Bottle (Anaerobic)		Urine Collection Bottle for Virus Isolation		
Blood Culture Bottle (Myco / F-Lytic)				
Blood Agar Plate		BIOCHEMISTRY SECTION		
Chocolate Agar Plate		Trace Element Serum Tube		
GC-lect Agar Plate (for GC Culture)		Trace Element K2 EDTA Tube		
Ureaplasma/Mycoplasma Culture Media		24hrs Urine Copper Container		
Sabouraud Dextrose Slant				
Enriched Thioglycollate Broth (TH)		CYTOLOGY SECTION		
Brain Heart Infusion Broth (BHIB)		Shandon Fixative (collection of		
		CSF/Vitreous fluid for Cytology		
		<u>Investigation only</u>)		
Transwab				
Teflon-Coated Slide				
Spore Strip (Bacillus atrophaeus)				
Spore Strip (Geobacillus stearothermophilus)				
Uricult Dipslides				
OC Light Kit for Stool Occult Blood				
TB Quantiferon Tubes (Four tubes per set)	(set)			

INSTRUCTIONS:

- 1) Kindly send request form one day in advance for the above order(s).

 Fax request forms (Fax no: 62254965) to Client & Specimen Management (CSM),

 Monday to Friday, 8.00 am to 4.00 pm
- 2) Requested items are to be collected the following day: 9.00 am to 12.00 pm.

For CSM staff use: Item/s issued by :	Date:
--	-------