

Make/Change Appointment Pay Online

## HUMAN HERPESVIRUS 6 (HHV6) QUANTITATIVE PCR

Synonym(s): HHV6 DNA Load

Back To Previous Page (javascript: history.go(-1))

Lab Section Category	Molecular Diagnostics - Microbiology  Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/moleculardiagnostics.aspx).
Indications	To detect, monitor and quantitate HHV6 DNA in plasma of patients following solid-organ transplants, bone marrow and stem cell transplants. In the absence of another identified etiology, detection of HHV-6 in CSF is also considered diagnostic of HHV-6 encephalitis.
Specimen Required	3-5 ml of whole blood in EDTA
Storage and Transport	Refrigerate sample until transfer to laboratory. <b>Do not freeze</b> . Send sample at 2-8 <sup>o</sup> C. Sample should reach the laboratory as soon as possible, preferably within 24 hours after collection.
Method	Real-time PCR (In-house)
Test Result	Detected (quantity reported in copies/mL), Undetectable or Inconclusive  Linear Range: 1,100 – 11,000,000,000 (3.04-10.04 log) copies/mL
Turnaround Time	5 days
Day(s) Test Set up	Tuesday and Friday

Back To Previous Page (javascript: history.go(-1))

**Change History Notes** 

28 Sep 2013 10:38 AM
 Day(s) Test Set-Up has been last updated on 28/09/13

 16 Jan 2017 11:10 AM
 Turnaround Time: Last updated on 16/01/17.
 ©Copyright 2018 Singapore General Hospital, Outram Road, Singapore 169608