SharePoint Adeline Sng Shi Hui v





## ALBUMIN, 24 HR URINE

Synonym(s): 24h Urine Albumin

Back To Previous Page (http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=Q2xpbmljYWwgQmlvY2hlbWlzdHJ5-kpqRmFbrvps%3d) (http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=Q2xpbmljYWwgQmlvY2hlbWlzdHJ5-kpqRmFbrvps%3d)

Lab Section Category	Clinical Biochemistry  Click here to find out more about the write-up (/clinical-departments-centers/pathology/pathology-handbook/lab-discipline-special-instructions/pages/clinical-biochemistry.aspx).
Specimen Required	24 hr urine (no preservative).  Urine preferably refrigerated during collection.
Method	Immunoassay
Reference Interval/Value	>12 yrs : < 30 mg/day (Normal) 30 - 300 mg/day (Moderate Albuminuria) > 300 mg/day (Severe Albuminuria)
Turnaround Time	Mean = 4 hrs, Range: 3 - 5 hrs
Day(s) Test Set up	Monday – Saturday

Remarks	Instructions before sampling:
	Patient should be well as baseline. They should have no urinary tract infection, no acute febrile illness, no
	intense exercise within the previous 24 hrs and not be menstruating.
	Total 24-hour urine volume must be clearly indicated.

Back To Previous Page (http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=Q2xpbmljYWwgQmlvY2hlbWlzdHJ5-kpqRmFbrvps%3d) (http://cms-sgh.sppub/Clinical-Departments-Centers/Pathology/Pathology-Handbook/Lab-Discipline-Special-Instructions/Tests/Pages/Forms/AllItems.aspx?category=Q2xpbmljYWwgQmlvY2hlbWlzdHJ5-kpqRmFbrvps%3d)

©Copyright 2019 Singapore General Hospital, Outram Road, Singapore 169608